



2017-2018 CATALOG

New & Refurbished Medical Equipment

Anesthesia Machines • Patient Monitoring
Surgical Lighting • Surgical Tables • Stretchers
Electrosurgical Units • Defibrillators/AEDs • Respiratory Ventilators
and much more



DRE Integra SP VSO2 • Page 5

- 4 Anesthesia
- 8 Patient Monitors
- 11 Central Monitoring
- 12 Exam Rooms
- 13 EKG Machines
- 14 Surgery Lights
- 18 Chairs & Tables
- 20 Surgery Tables
- 22 Mobile Surgery Tables
- 23 Stretchers
- 24 Infant Care
- 25 Fetal & Maternal Monitors
- 25 Ultrasound
- 26 Laparoscopy
- 26 Video Laryngoscopes
- 27 Microscopes
- 28 Defibrillators & AEDs
- 30 Electrosurgical Units
- 32 Ventilators
- 34 Syringe Pumps
- 35 Infusion Pumps
- 35 Aspirators
- 36 Warmers
- 37 Hospital Beds
- 37 Autoclaves
- 37 Stainless Steel
- 38 Service & Repair
- 39 Supplies and Accessories
- 39 Medical Gas



Introducing DRE's

Outfit your operating room quickly and economically, with peace of mind that your investment will last.

Centrally source all the the O.R. equipment you need with Turn-Key Equipment Solutions from DRE. You can trust our expert sales representatives to provide all the new and refurbished medical equipment you need for the lowest possible price.

PAGE 26 Stryker Laparoscopy Tower

Visualization and Documentation System



PAGE 36 DRE Comfort Series

Blanket and Fluid Warmers



PAGE 14 DRE Maxx Luxx 160 Surgery Lights



PAGE 4 DRE Integra SL Anesthesia System



PAGE 9 DRE Waveline Pro

O.R. Monitor available with EtCO₂ and 5-agent analysis



PAGE 28 Physio-Control Lifepak 20

Defibrillator and Monitor



PAGE 35 DRE DM-660 Aspirator



PAGE 22 DRE Versailles P100 Mobile Surgery Table



Operating within your budget for more than 30 years

CALL TODAY FOR A QUOTE!

1-800-499-7082 | dremed.com/2017

Every day, healthcare professionals rely on medical equipment from DRE to provide quality care to their patients. Your facility can save time and money by turning to DRE for all the new and refurbished equipment you need. DRE is:

A nationally-recognized medical equipment partner for more than 30 years.

Since we began in 1984, our company has established itself as a knowledgeable equipment provider with extensive industry knowledge. You can trust us to know the market, and to guide you to the right products for your practice.

A single supplier for all of your medical equipment needs.

Trust our friendly staff to provide affordable and innovative solutions for your practice. We strive to be a complete equipment resource, providing each customer with a pleasant equipment procurement experience.

Two Cost-Saving Medical Equipment Solutions

1 New DRE-brand Equipment

Our comprehensive lines of medical equipment are dependable and include the innovative features needed in the modern medical environment — at economical prices that fit nearly any budget. We offer a variety of warranty options, protecting your investment and providing you with peace of mind.

Advanced & Economical Factory New Equipment



2 Certified Refurbished Equipment

For more than 30 years, DRE has been a leading refurbisher of medical equipment. Our wide selection of refurbished equipment includes patient monitors, anesthesia machines, defibrillators, electrosurgical units and surgical tables — including current-model refurbished equipment with advanced features.



DRE Integra SL

Electronic flow display with digital and traditional flowmeters.

Fully integrated absorber ventilator with 6 ventilation modes.

Illuminated work space, large capacity drawer units and top shelf with secure GCX mounting system for patient monitors.

Fully compliant to ISO 80601-2-13 and RoHS standards.



< Patient monitor sold separately

Intuitive design in a compact footprint.

Web ID 2350

FIND IT FAST at www.dremed.com when you search by Web ID number

Patient monitor > sold separately



DRE Integra SL3

Selectatec-compatible backbar can fit up to three vaporizers.

Electronic flow display and conventional flowmeter options.

Fully integrated absorber and ventilator.

Hospital Information System (HIS) capabilities when using electronic flow display or VFDBullet.

Expertly designed, clinician-focused machine for high traffic operating rooms.

Web ID 2403

DRE Ventura CS5

High-quality anesthesia workstation with a 10.4" TFT color touchscreen with large, easy-to-read controls.

Fully-integrated, ergonomically positioned ventilator and absorber.

Multifunction anesthesia ventilator designed for adult and pediatric use. Advanced ventilation modes include volume-control and pressure-control.

A safe, reliable anesthesia system built for the modern OR.

Web ID 12520



DRE Integra SP VSO2

Small compact unit with retractable top handle for easy transport.

Three gases (O₂, N₂O and air).

Two-year manufacturer's warranty.

Mechanical Anti-Hypoxic Device and Air/N₂O interlock.

Conveniently located side gas ports.

Selectatec®-compatible vaporizer mounting system.

Compact system is great for office-based anesthesia.

Web ID 40



Available with an optional absorber

Available with a hard-shell carrying case, making the system easier to transport.

Available with a durable mobile cart.

New and Recalibrated Vaporizers

Penlon Sigma Delta (new) delivers accurate concentrations under many conditions, even at low flow.

Recalibrated options include Datex-Ohmeda Tec Series, Dräger Vapor 2000, 19.1 and 19.3.

Clean and Calibrate your vaporizers with DRE!

TEC 3, TEC 4, TEC 5, Dräger 19.1, Ohio and Penlon.

Also have your TEC 3 and TEC4 vaporizers converted from Halothane to Isoflurane.

Penlon Sigma Delta

Web ID 56



DRE ANESTHESIA PARTS AND SERVICE



Certified Biomedical Equipment Technicians
Trained, Educated, Experienced

Factory-trained with Certification to Service Anesthesia Systems from GE/ Datex-Ohmeda, Penlon, Dräger, DRE, and more.





Dräger Apollo

E-Vent plus servo controlled high-speed piston ventilator provides performance previously only seen in the ICU.

Designed to support the user in low and minimal flow settings.

Simple, straightforward start-up reduces the burden on the staff and provides complete information on device status

Easy-to-use anesthesia system with user-friendly ergonomics.

Web ID 2120

Open-Platform Architecture
can easily integrate into overall workflow



Dräger Fabius GS

High-contrast color monitor with Dräger user interface for easy and familiar operation.

All major ventilation modes.

New, highly maneuverable trolley with central brake.

LED workplace illumination, spacious design, and large drawers.

Advanced ventilation technology with reliable and proven breathing system components.

Web ID 2248



Dräger Fabius Tiro

Intelligent safety features for enhance patient protection.

Complete with advanced features such as fresh gas decoupling, which maintains constant tidal volumes regardless of fresh gas flow changes.

High-quality ventilation in small spaces.

Web ID 2213

Open-Platform Architecture
can easily integrate into overall workflow



Datex-Ohmeda Aisys

Ventilation, vaporization and gas delivery are electronically controlled and monitored.

Data capabilities enable digital interfacing of patient and Carestation information.

Easily upgradable to enhance the functionality, add support, and enhance data management.

Provides digital feedback on ventilation, drug levels, patient status, and fresh gas usage.

Electronic protocols provide data continuity.

Precise digital controls ideal for low flow anesthesia.

Web ID 2204



Datex-Ohmeda Avance

Compact anesthesia system features advanced anesthesia delivery and premium anesthesia patient monitoring.

Clinical information management systems, that can include DEIO record keeper or Centricity Anesthesia.

Web ID 2203

Datex-Ohmeda Aestiva/5

Optimized for low flow anesthesia and reduced agent consumption.

Open systems architecture and lower overall height.

Web ID 1469



Datex-Ohmeda S/5 Aespire

Enhanced integration for Datex-Ohmeda monitors.

Lightweight and compact for seamless maneuvering.

Optional integrated auxiliary O₂ flowmeter and suction control.

Available with the advanced 7900 SmartVent — includes sophisticated ventilation capabilities, helping you meet the needs of neonates, trauma, compromised and routine patients.

Modern economical anesthesia.

Web ID 1818



DRE Waveline Touch

Bright, easy to use, touchscreen display with up to 8 adjustable waveforms.

Optional capnography and anesthesia gas.

Wireless and central station compatible. USB and network access to software and data.

Multi-lead ECG monitoring with optional EtCO₂.

Web ID 1979

Parameters: ECG, Pulse Oximetry, Noninvasive Blood Pressure, Heart Rate, Temperature, Respiration.

Optional: EtCO₂, Dual Invasive Blood Pressures, Anesthesia Gas, 5-agent Analysis, and FiO₂.



Vibrant 12.1" user-friendly touchscreen for seamless operation

DRE Waveline EZ

Transportable monitor with an 8" touchscreen.

Long-lasting battery back-up (up to 4 hours).

Optional CO₂ monitoring.

Capnography available.

Web ID 1980

Parameters: ECG, Pulse Oximetry, Noninvasive Blood Pressure, Heart Rate, Temperature, Respiration.

Optional: EtCO₂ and Printer



Compact design convenient for mobile monitoring



DRE Waveline EZ MAX

10.4" Color TFT LCD touchscreen monitor.

Arrhythmia, S-T Analysis, OxyCRG, large font display.

Non-proprietary CO₂ technology saves you consumable expenses.

Web ID 2402

Parameters: ECG, Pulse Oximetry, Noninvasive Blood Pressure, Heart Rate, Temperature, Respiration.

Optional: Nellcor OxiMax™ SpO₂, 2-IBP, ICP, C.O. Respiration CO₂ (Mainstream and Sidestream)

DRE Waveline Pro

Available with a built-in anesthetic agent bench!

Large (15"), touchscreen with high-resolution flat panel display.

Optional Integrated Dräger Anesthesia Gas Bench expands measurement capabilities.

Wireless and central station compatible. USB and network access to software and data.

Displays as many as eight waveforms.

Touch Screen O.R. monitor available with EtCO₂, anesthesia gas and 5-agent analysis.

Web ID 1821

Parameters: 3-lead or 5-lead ECG with arrhythmia detection, Noninvasive Blood Pressure, motion tolerant SpO₂, Temperature, and Respiration rate.

Optional: EtCO₂, Dual Invasive Blood Pressures, Anesthesia Gas, 5-agent Analysis, 5-agent ID and FiO₂.



GE Dash

Complete Series of Models Available



Web ID 1952

GE CARESCAPE B850



Web ID 12527

GE Datex-Ohmeda CardiCap/5



Web ID 1691



Philips IntelliVue Systems

Complete Series of Models Available

Monitors and central station components.

IntelliVue MP90, 70, 50, 40, 30, 20, & 5 Telemon Units Also Available

Web ID 2158

High-performance monitoring for flexible care settings.





DRE Echo Series

Cost-effective vital signs monitor.

5.7" high resolution display clearly presents vital sign information.

Optional thermometer for obtaining temperature in approximately 15 seconds.

Available with your choice of Pulse Oximetry, Noninvasive Blood Pressure and Covidien Quick Temp. Wireless Capable.

Ultra-lightweight and portable with a 10-hour battery.

Web ID 2148

Suitable for adult, pediatric and neonatal patients

Welch Allyn VSM 300

Affordable vital signs monitor that is both flexible and easy to use.

Settings for adult, pediatric and neonatal support.

Wireless Capable.

Web ID 449

Other Models Available



ProCare 400

Easy-to-use features and an intuitive menu structure.

NIBP, SpO₂, Temperature, and a Recorder.

Wireless Capable.

Web ID 1726

Other Models Available



Certified Refurbished

CO₂ Monitor • Hand-held • Bedside

DRE Echo CO₂

Respironics Loflo sidestream or CAPNOSTAT 5 EtCO₂ mainstream measurement.

Real-time parameters measurement display with trend table for easy reviewing.

Bi-directional communications with central station by wired or wireless network.

Built-in rechargeable Lithium-ion battery for 10 hours continuous operation.

Lightweight, portable design and user-friendly interface.

Flexible configurations to meet different clinical needs.

Standard SpO₂.

PR measurement.

5.7" high resolution display for easy reading

Effective capnography monitoring.


Web ID 2150



DRE Envoy TS

Monitor up to 16 patients throughout your facility from a single location with the DRE Envoy TS Telemetry System.

Web ID 1235



ECG Transmitters

Waveline EZ

Waveline Touch

Waveline Pro

The Envoy TS integrates with ECG transmitters and DRE patient monitors

The Envoy TS is economical and configurable, and features telemetry transmitters for wireless cardiac ECG monitoring.

No other system provides the high level of flexibility offered by the DRE Envoy TS.


Each channel is independent, allowing the use of stationary, wireless transport and ECG telemetry transmitters on a single system.

Operating system is Windows-based and provides a highly versatile platform.

Supports network interfaces and laser printers.

HL7 Output for EMR Compatibility.

Monitor Patients from Pre-op Through Post-op



Easy-to-use interface

Split-screen technology allows you to review, edit and enter patient information while observing all patients in the system.



ECG Transmitters **Small, portable and lightweight**

Single, multi-lead

Frequency: Medical 608 MHz

Independent lead wires

Waterproof

Operates on two AAA batteries

Flexible and durable



ONE FURNISHED ROOM • ONE AFFORDABLE PRICE • ONE GREAT DEAL

DRE ONE Room

Practice Exam Room Package

Four great products, packaged at one affordable price.


Furnish an entire exam room with contemporary style in one convenient package. The One Room from DRE provides a modern base and wall cabinet combination, a versatile exam table, comfortable pneumatic stool and side chair, all color coordinated and ready to go.

There is no better exam room package of equipment at a better price point than The One Room. The equipment alone, backed by a **5-year warranty**, is a big deal. The affordable price makes it a great one.

Web ID 2398

Colors available:
Sapphire, Cream and Black
Made in the USA

Optional Equipment:



DRE Vista LED
Examination Light



DRE Echo CO₂
Vital Signs Monitor



DRE WS350
Wall Mount Diagnostic System



DRE Patient Exam Table

Versatile and durable with multi-use features. Perfect for the private practice physician or clinic.



DRE Cabinet Combo

Features stainless steel sink, all laminate construction, soft close hinges and drawers that hold up to 75 lbs.



DRE 5-Leg Pneumatic Stool

Easily adjustable; 15" diameter, padded seat and dual wheel casters for durable performance.



DRE Premium Side Chair

Features a durable, metal, powder-coated frame and comfortable padded seat and backrest.

DRE True ECG Advance

12.1" color touchscreen.

Real-time waveform freezing.

Automatic measurement and interpretation tested with CSE database.

Alphanumeric keyboard and one-touch operation.

200 ECGs in internal memory.

Review 120 seconds of ECG waveform.

Touchscreen 12-channel ECG.

Web ID 1977



Available with optional cart

Optional stress testing can support treadmills and bicycles

DRE True 12T

8.4" color touch-screen.

Preview report for both real-time EKG printing and file storage.

Optional EKG printing styles, printing speed, sensitivity, and filtered waveforms for less repeat examinations.

Built-in rechargeable Li-ion battery.

Excellent cost performance with large touch screen display.

Web ID 1987



Compact, lightweight design with carrying handle

Philips TC70

15" color touchscreen (fast and easy).

Philips DXL 16-Lead ECG Algorithm provides industry standard diagnostic capabilities for a wide patient population.

State-of-the-art performance.

Web ID 2440



DRE Velocity Pro EKG

Comprehensive software with multiple applications.

Multi-user permission configuration.

Auto arrhythmia detection.

Optional stress test function.

30s-360s HRV analysis.

Compatible with Windows 8 and 10.

PC-based 12-channel ECG system.

Web ID 1947



INSTALL
the software
CONNECT
the Velocity EKG
MEASURE
ECGs on your PC

Visit dremed.com for complete hardware requirements.





DRE SLS 9000 Plus

Maximum output rated at 140,000 Lux.
 Rated for 50,000 hours of LED life.
 Expertly engineered design makes it easy to position and adjust light heads for any surgical procedure.
 LCD touch screen operation measures illumination intensity, color temperature and Chromogenic index.

Bright, consistent illumination.

Web ID 12686

Available as a single or dual head ceiling mount

DRE Vision LED

Produces unsurpassed illumination with economical energy consumption
 50,000-hour rated LED life.
 100,000 lux output at 1 meter.
 4,300° Kelvin color temperature.
 5-stage dimming and on/off switch on the patented sterilizable handle.

Provides unsurpassed longevity and flexibility.

Web ID 2119



Efficient LED lights provide low heat

DRE Maxx Luxx LED 160

LCD touch-screen for easily adjustable output.
 Cool output LED light source.
 Built to last at a rated 50,000 hours.
 160,000 lux output (14,800 foot candles).
 Adjustable focusing and CRI.
 Adjustable color temperature for optimizing human tissue colors.
 Available in dual ceiling, single, and wall configuration
 Web ID 2393



Available in multiple mounting configurations including single and dual wall or ceiling mount and a portable floor model



DRE Maxx Luxx LED

Highly efficient surgery lighting system.
 Extremely cool output — lowest possible heat.
 Superb color temperature of 4,300° Kelvin.
 Extremely long life at a rated 50,000 hours.

Maximum brightness with extremely cool output for 50,000 hours.

Web ID 2295

Available in multiple mounting configurations with light heads, camera, and monitor combinations
 Portable stands also available with a light head



DRE Vision EX5

Extremely cool output — lowest possible heat.
 75,000 lux output at 1 meter.
 Extremely long life at a rated 50,000 hours.
 4,300° Kelvin color temperature.
 LED array delivers maximum sharpness and color accuracy.
 3-stage dimming.

Midsized LED surgery light offers efficiency and economy.

Web ID 2400

Available in multiple mounting configurations including single and dual wall or ceiling mount and a portable floor model





DRE Vista LED

Maximum output rated at 50,000 lux.
4,300° K bright white illumination.
Life rating of 50,000 hours.
2-stage dimming.

Compact head design with smooth arm articulation.

Offers brilliant illumination at an affordable price.

Web ID 2196

Provides low, economical energy consumption

Available in dual ceiling, single, wall and mobile floor model configurations



DRE FX-300+ Xenon

100% brightness at startup. No lag or warm up needed.
1,000-hour expected xenon lamp life is one of the industry best.
The rotating turret accepts cables made by ACMI, Wolf, Storz, Olympus and more.

Web ID 2128



High performance 300-watt xenon light source with optional video system.



DRE VideoFX

In Procedure Video Capture

Through the Lens Headlight Video System captures procedures from the surgeon's line of sight

Record and play back cases with standard composition or S-Video output to your preferred recording device.

Web ID 2222



DRE FX-10 LED

10-watt LED light source with a long life: 25,000 hours.
Portable, light and compact; weighs less than 10 ounces

Versatile: Single and Tri-Port models available.

Firmly mounts to steel surfaces magnetically.

Web ID 2219



Cool-running LED light that starts to full power instantly.

Surgical Headlights

DRE Xavier-C10 LED

Operating room quality, high intensity 10-watt LED Headlight rated for 25,000 - 50,000 hours use.

Rechargeable battery pack guarantees 4.5 hours of maximum intensity light.

89,000 lux at distance of 30cm, 55,000 lux at distance of 50cm.

Clear field of illumination prevents hot-spots and shadows.

Web ID 2380



DRE Xavier C3 LED


Portable light-weight LED.

Rechargeable battery; 7.5 hours of high intensity light.

49,000 lux at 10" working distance.

Web ID 2293

Halogen model upgrade available!



DRE Xavier 3 LED

Light-weight; cable free.

80,000 lux at 15 cm.

Easily recharge Li-ion battery in 2.5 hours.

Web ID 1943



DRE Vista Flex LED

Cost-saving, energy-efficient LED bulb produces up to 50,000 hours of illumination.

Provides a color temperature of 5,300 K.

Light produces 70,000 Lux at an eight inch working distance.

Control panel features intuitive controls and two handles for easy maneuverability.

Available in a variety of mounting options including mobile stand, wall mount or table mount.

Versatile light with maximum positioning capabilities.

Web ID 12689



DRE SLS 2500 LED

Shadowless lamps provide color temperature close to natural sunlight with 3,500 ~ 5,000 K.

60,000 lux output at 1 meter.

Light and strong nova alloy arm allows easy and convenient positioning.

Web ID 2351



DRE Vista PRO 3 LED

Alternate between 3 preset color temperatures: 4,400°, 3,800°, and 3,000° Kelvin.

40,000 lux output at 0.5 meters.

Web ID 2407



Sleek Chair with a Small Footprint




DRE Milano E20

Power procedure chair naturally adapts to almost every physician, patient, procedure and environment.

Comfortable for both the operator and patient: the Milano E20 is wider than traditional procedure chairs and lowers to 20.5", helping patients with mobility difficulties get on and off the chair.

Foot pedal and handset remote adjusts the height, back, Trendelenburg and level positioning with 4 actuators and includes an automatic reset control.

Web ID 2199



Multiple color options available



DRE Milano T50

A universal power procedure table with fully-electric multifunctionality. Includes four actuators and standard foot pedal control.

The electronic control box and the actuators are versatile and easy to service. Featuring an automatic reset switch that returns the table to its starting position.

Available in 110 and 220 V cc. configurations.

The galvanized steel frame is fully covered with ABS for an attractive appearance as well as safety.

Designed for Dermatology, Oral, Plastic & Cosmetic Surgery, OB/GYN, ENT, and General Procedures.

Web ID 2200



DRE Milano OB20

Full electric operation OB/GYN Procedure Chair.

Simple, automatic operation with 3 independent movements: height, backrest and Trendelenburg.

20" minimum height ensures easy access for the patient; 40" maximum height ensures optimum position for the specialist.

90° rotation to both sides and adjustable headrest.

Web ID 2339



DRE Milano D20

Wheeled base for moving the chair throughout your facility.

Articulating headrest can be adjusted for a variety of procedures.

Ideal for dental procedures, implantology, orthodontics, ENT, facial procedures, examinations, diagnostics and minor surgeries.

Four Independent Electronic Movements: Height, Backrest, Leg Rest and Seat Rest.

A highly-adjustable, affordable chair for modern practices.

Web ID 12499



Oakworks DTPM300



3 movement imaging table solution.

Trendelenburg/Reverse Trendelenburg function, complete C-Arm compatibility and locking casters.

Web ID 1241

Oakworks CFPM400



4-movement C-Arm Imaging Table.

Low height for patient transfers to maximum height. C-Arm mobility requirements for full lateral imaging.

Web ID 1663

DRE Milano P20

Feature rich, stylish and affordable Podiatry Chair.

Leg supports swivel independently up to 50°.

Foot pedal adjusts the 3 electric movements: height, backrest and seat tilt.

30° Partial Trendelenburg positioning.

Web ID 2395



DRE Milano H50



Ideal for: cosmetics, hair transplantation, ophthalmology, any lengthy outpatient procedure.

Hours of breathable comfort for long procedures.

Web ID 2415

Midmark 222 & 223

Barrier-free high-low power exam table.

Low and high table positioning.

Web ID 1154



Operating within your budget!



DRE Torino 550

Multi-purpose positioning capabilities.

Longitudinal slide up to 13.78" for easy C-arm access.

Table weight load capacity rated at 550 lbs.

Adjustable height from 28" to 40" for good access to the operative field.

Table movements include: Trendelenburg, reverse Trendelenburg, leg and back plate adjustment and lateral tilt.

Large casters and convenient braking system for easy and convenient maneuvering throughout any surgical facility.

Battery backup for continued use in the event of a power loss.

Web ID 12645



Memory padding is waterproof and antistatic



Optional Accessories:
Kidney bridge
Shoulder support
Arm support
Long fixing clamps
Leg and foot supports



DRE Torino EXL

Features long horizontal sliding, C-arm and x-ray access.

Micro touch remote control with low noise movements.

Flexible adjustments on head plate, back plate, seat plate and built-in kidney bridge.

Quiet and reliable electric motors.

Stainless steel design for excellent stability.

550 lbs weight capacity, 25° Trendelenburg, 15° Lateral Tilt (left, right).

Web ID 2254



Excellent choice for a variety of operations



Included Accessories:
Anesthesia screen
Shoulder support
Body support
Arm support
Leg support
Foot plate

Steris 4085

Four-section tabletop contours to virtually every patient posture.

Removable leg section provides unrestricted perineal access.

Powered, radiolucent kidney elevator for lateral procedures.

Generous height and weight ranges accommodate many patient types.

Compensates for uneven OR floors with unique self-leveling floor locks

Provides outstanding access for C-arm.

Web ID 2130



Certified Refurbished

Steris Amsco 3085 SP

Safely raise and lower larger patients in both normal and reverse orientations.

Unprecedented versatility for all surgical specialties.

Designed for flexibility in patient positioning.

Exceptional image amplification (IA) coverage.

Durable, flexible surgery tables support up to 1,000 lbs.

Web ID 1499



Certified Refurbished

Skytron 6600

General use electro-hydraulic surgical table.

Can lift and fully articulate patients weighing up to 1000 lbs.

A wide selection of models available.

Web ID 1755



Certified Refurbished

Skytron 3500

21" top slide, 500 lb lift/articulation, return-to-level and return-to-center.

Provides full body imaging.

Satisfy all of your C-Arm and general surgery needs.

Web ID 1754



DRE Versailles P100

The Versailles P100 is a complete surgery solution that can move patients from the pre-operative area to the operating room and then to the recovery area.

The Versailles P100 does it all by incorporating stretcher mobility into a multi-purpose surgery table!

Lightweight and maneuverable.

75% C-arm access.

Increased surgical access.

Low table height — helps patients more easily get on and off the table.


Manual option available - DRE Versailles M75.

Web ID 2183

Carries patients from pre-op to the O.R. to recovery — no patient transfers required!



Surgery Table Positioning and the Mobility of a Stretcher



Ophthalmology

Patient platform traversed to head end for access.



Gynecology

Patient platform traversed to foot end for access.



Arthroscopy

Optional arthroscopy backrest fitted.



C-arm Coverage

Powered configuration provides full C-arm coverage.



Headrests available for surgical procedures

DRE Millennium 5

Backrest lifting controlled by oil hydraulics.

Knee-rest lifting controlled by adjustable crank.

Central braking system with steering function.

Steel bed frame with powder coating.

3" high density PU foam mattress.

Collapsible side rails on both sides that fold manually.

8" steering wheel central brake on four casters.

Web ID 2089

Durable patient stretcher available with a fifth wheel system for improved steering.



Two integrated I.V. poles



Designed for high volume hospitals and surgery centers

DRE Savoy QA3

Designed for maximum ease of use combined with superior patient comfort.

Vertically operated patient safety rails that stow away under mattress for a virtually zero transfer gap.

Fixed transfusion pole with quick release locking mechanism.

Hydraulic lift with low height capability adjusted by pedals located on both sides of the stretcher.

5th wheel steering for easier maneuverability.

Gas assisted backlift for enhanced user and patient comfort.

Trendelenburg tilt allows rapid tilt to be achieved even at lowest height.

K8 pressure care mattress manufactured from temperature sensitive gel.

Web ID 2270

Ergonomically positioned pushing handles



Smooth and controlled Trendelenburg tilt





GE Giraffe OmniBed

Protects from heat loss while the bi-directional airflow maintains thermal stability.
 Infant receives uniform heat regardless of position.
 Radiant warmer provides consistent heat in any position.
 No need to transfer patient from incubator to infant warmer.
 360° rotating patient mattress.

Full-featured incubator and a radiant warmer in one.

Web ID 2330



Air-Shields C-2000

Bi-directional air flow for consistent heating and humidifying.

Optional Dual Access Doors – Allow access from either side, expediting care.

Large mattress and dual probes allow for simultaneous care of twins in the same unit.

Web ID 1738



Air-Shields TI-500

Mobile neonatal intensive care system designed for transport.

Front mounted controls for quick access.

Auxiliary flow meter with an additional hand circuit to allow manual resuscitation.

Web ID 1739



Dräger Babytherm

smartSWIVEL mechanism focuses radiant warmth on infant.

Web ID 1735



Dräger Resuscitaire

Precise gas delivery prevents respiratory complications in neonates.

Web ID 2202



Integrates with DRE-brand Fetal & Maternal Monitors

DRE Envoy FM

Central station for fetal/maternal monitors.
 Provides continuous trace values and displays all measurement values for the fetus and mother.
 Provides audio and visual alerts, making doctors and nurses immediately aware of critical events.
 Saves time and effort with a user-friendly interface.

Central station manages up to 32 fetal/maternal monitors at once.

Web ID 1995

DRE Compact FM



Dual fetal heart rate and TOCO monitoring.
 10.2" folding color wide-screen.
 Multi-crystal pulsed wave transducer.

Web ID 1940

DRE Compact II



5.5" folding graphics display screen with large, easy-to-read numerals and waveforms
 Backup memory for seamless monitoring.

Web ID 2394

GE Corometrics 250cx

Comprehensive perinatal monitoring needed during labor and delivery.

ChartLight with backlit keys and dimmable display

Web ID 12566



GE Corometrics 171/172

Monitors external uterine activity

Highly portable unit at just 8lbs.

Dual ultrasound capabilities for twins pregnancies

Web ID 2459



Ultrasound Systems

GE Logiq Book

Full ultrasound capabilities in a 10-pound, hand-carried, battery-powered scanner.

Web ID 1825



DRE Crystal 4PX

Next generation imaging and software.
 High resolution color LCD screen.
 Lightweight mobile design.

Easy-to-use patient management system for storing DICOM files.

Ultra-sensitive Spectral Doppler with High Pulse Repeat Frequency (HPRF).

Web ID 2408



Laparoscopy Towers



Endoscopy Digital Camera
Stryker 1188 HD and 1288 HD available.

Light Source
Stryker X-8000 and L-9000.

Digital Capture System
Stryker SDC Pro 2 DVD.

Insufflator
Stryker 30L, 40L, and pneumosure.

Web ID 1866

Light Source
Olympus CLV-190, CLV-180, and CLV-160 available.

Video System Center
Olympus CV-160, CV-190, and CV-180.

Insufflator
Olympus UHI-2, UHI-3, and UHI-4.

Web ID 12718



Certified Refurbished

Stryker Endoscopy 1288 HD

3 CHIP digital camera
Also available: Stryker 1088 HD and 1188 HD



Olympus CV-180

Equipped with high-resolution HDTV and Narrow Band Imaging (NBI)



Olympus GIF-XP160 Video Gastroscope
5.9 mm insertion tube



Olympus GIF-N180 Video Gastroscope
Large, sharp images

Video Laryngoscope

DRE InView HD

Video laryngoscope ideal for a wide range of applications: ICU, Cardiology, Emergency and Ambulance.

Lightweight and ergonomic design weighs less than one pound.

3.5" display provides high-performance visualization capabilities for every intubation.

Low-cost disposable blades gives you huge savings.

Battery life lasts approximately 200 minutes.

Video display rotates 180 degrees, as well as up and down.

Portable and ergonomically designed.

Web ID 12688



Designed to provide an excellent view of the glottic opening and improve the success of tracheal intubation.



Zeiss Opmi Pico ENT

Zeiss optics: Outstanding focusing abilities with razor-sharp images; Homogeneously illuminated field of view; 5-step, apochromatic magnification changer; Focusing objectives: f = 200 mm, f = 250 mm; 200 mm to 400 mm in steps of 50 mm; 12 V 100 W illumination.

Web ID 12714



Certified Refurbished

Zeiss Opmi MD

Zeiss optics: OPMI MD Microscope Body, f=175 mm Objective Lens, 0-60 Degree Inclinable Binocular, 150w Halogen Illumination, Foot Switch Control (Zoom and Focus).

Web ID 1808



Leica M-690

Leica optics: Excellent Depth of Focus, Red Reflex Enhancing Objective Lens, High Resolution, Internal, Motorized Focus and Zoom, 0-180 Degree Inclinable Binocular, 1:5 Magnification Ratio.

Web ID 332



DRE Compass LED

Surgical Microscope

Combination of classic aesthetics and modern technology



Web ID 2290


MOUNTING OPTIONS:

- Ceiling
- Wall
- Fixed stand
- Table top clamp
- Floor mount



0-180° Inclinable Binocular

This ergonomic feature allows for comfortable and exact microscope positioning.



High Quality Lens

German Optical Glass, high resolution, real 3-dimensional precise image reproduction, large depth of field, wide field of view, and good contrast



6-Step Magnification

provides a larger magnifying range and more selection on magnification to meet the requirement of your practice.



Integrated LED light source

provides long-lasting, pure white color light.

Physio-Control Lifepak 20

Lifepak 20 has seven main operating modes.
 AC line operated device with an internal battery as backup.
 Fully featured defibrillator/monitor (pacing and SpO₂) 5.58 kg (12.3 lbs).
 Displays a minimum of 4 seconds of ECG and alphanumeric for values, device instructions or prompts. Option to display one additional waveform.
 Web ID 1681



We test, touch up, replace, and restore these units to manufacturer specifications and offer a great value on recent technology.



Philips HeartStart XL

A combination of manual defibrillation and automated external defibrillator capabilities.
 Enables the first caregiver on the scene, whether an ACLS or BLS clinician, to deliver lifesaving defibrillation therapy.
 Web ID 1487



SMART Biphasic Technology

Clinically proven superior to high-energy, monophasic waveforms for efficacy and for minimizing post-defibrillation heart dysfunction.

Quick Charging

Charges to highest energy level, 200 Joules, in less than 3 seconds.



Zoll M Series

Lightweight for easy handling — the built-in AC power/battery charger does not add size or weight.
 Zoll Uniform Operating System features consistent, straight-forward controls; and intuitive menus lessen training costs, minimize operator confusion, and reduce the chance of error.
 Zoll's proprietary Rectilinear Biphasic low-energy waveform.
 Web ID 1293

Built with Patients and Users in Mind



Physio-Control Lifepak 12

Component-based and upgradable system is designed for long-term use.
 Standard paddles, pediatric paddles, and/or sterilizable internal and external paddles available.
Compact and portable defibrillator/monitor — weighs only 14.5 lbs.
 Web ID 1435



Philips HeartStart MRx

Superior diagnostic measurements and a patented SMART Biphasic resuscitation waveform in a single, thoughtfully designed, lightweight device.
 Monitoring and therapy data are clearly and logically arranged on-screen.
 Adjustable ECG size and autogain



ALS monitor/defibrillator with CPR measurement and feedback technology.

Web ID 1492



Zoll E Series

Patient Connection: 3-lead, 5-lead or 12-lead ECG cable, paddles or MFE pads, selectable by front panel switch.



Designed for EMS use

Web ID 2343

Zoll R Series

Pacing, end-tidal CO₂, non-invasive blood pressure, and SpO₂ monitoring capabilities.
 Web ID 2383



Cost-effective advanced resuscitation

Zoll AED Plus

Gives even infrequent rescuers the confidence they need to help save lives.
 Simple audio prompts and pictures reinforce every step in the resuscitation process.
 Real CPR Help with instant feedback on rate and depth of compressions.
 Easy-to-use CPR-D-padz for fast, accurate placement.
 CPR-D-padz is based on extensive human anthropometric data and studies designed to accommodate the wide range of patient sizes and shapes.



Step-by-step support for the entire chain of survival.

Web ID 1182

One of the best values in the market. These dependable generators offer advanced technology, smooth cut and a wide range of features.



Optional durable cart available



DRE ASG-300

Specifications (in watts)
Cut: 300 W | Blend: 200 W
Pinpoint: 120 W | Spray: 80 W
Bipolar: 80 W

Delivers consistent, repeatable power into varying load impedances

10 blend setting

Defibrillator-proof type CF equipment

Web ID 244

Split or solid return electrodes • No calibration required • One-year warranty available

Upgrade available with dual monopolar and tissue fusion technology



Bovie Smoke Shark II

35-hour filter combines filtration in one casing, keeping cost down and minimizing waste.

Filter captures particles down to 0.01m with 99.99995% efficiency.

24 ounces of carbon for superior odor elimination.

Quiet, effective smoke evacuation.

Web ID 2409

Bovie Derm 101

10-watt high frequency monopolar desiccator to perform in-office, minor skin-related procedures.

Bright LED display for clear readings anywhere in the procedure suite.

Web ID 2406



Bovie Derm 942

Digital desiccator with 40 watts of power.

Power dial is adjustable within 1/10th watt increments.

Web ID 12693



Valleylab ForceTriad

Provides the operating room with one energy system for electrosurgical and tissue fusion needs and gives surgeons the ability to deliver precisely controlled energy and improved tissue effect.

First-of-its-kind integrated solution for open and laparoscopic electrosurgical procedures.

Web ID 2139



Certified Refurbished



Valleylab Force FX-C

For monopolar and bipolar procedures.

Power for cutting, desiccating, and fulgurating tissues.

A smoother cut through all tissue types.

Web ID 1436



Optional durable cart available



Valleylab Force 2

Blended cut modes provide flexibility through varying degrees of hemostasis.

300-watt output in cut setting.

Web ID 107



Certified Refurbished

Valleylab Ligasure

Revolutionary technology replaces almost all other hemostatic tools — it actually fuses vessel walls to create a permanent seal.

Web ID 1683





Rugged & Reliable for a Low Cost of Ownership

Puritan Bennett 840

High-performance pneumatics, dual-microprocessor electronics and DualView™ touchscreens.

Supports communication with patient monitoring and information systems, and provides seamless electronic data transfer to a patient's medical record.

Enhanced noninvasive ventilation (NIV): improves comfort for patients and performance for clinicians.

Delivers sensitive, precise breaths to critically ill neonatal, pediatric and adult patients.

Web ID 1236



Certified Refurbished

Puritan Bennett 980

Engineered to help patients breathe more naturally.

Versatile respiratory ventilator offers non-invasive SIMV and CPAP.

Tube compensation accurately overcomes the work of breathing imposed by the artificial airway.

Standby mode pauses ventilation and maintains settings while the patient is disconnected and auto-detects reconnection.

Networking capabilities allow interfacing with existing EMR systems.

Web ID 12683



Siemens Servo-i

Flexible ventilator is adaptable for all patient categories.

Four configurations are available: SERVO-i Infant, SERVO-i Adult and two editions of SERVO-i Universal.

SERVO-i supports both invasive and non-invasive ventilation.

Ease of training, operation and maintenance as well as flexibility and modularity ensures a low life cycle cost.

Web ID 1674

High performance, mobile, cost effective respiratory system.



Certified Refurbished

Respironics Esprit

Microprocessor-driven unit offers a range of modes and breath types.

Mobility and versatility make the Esprit ideal for use in any setting.

Web ID 2117



Pulmonetic LTV Series

Ventilator pre-sets allow for the quick set-up of new patients.

Automatically configures the initial ventilation settings for infant, pediatric or adult patients.

Intuitive interface features color-coded controls and alarms.

Invasive and Non-Invasive modes of ventilation for patients as small as 5kg.

Features an innovative Spontaneous Breathing Trial.

Web ID 1875

Ease of use and versatility help ensure optimal patient care in clinical settings.

Visit dremed.com for more LTV Series ventilators



Certified Refurbished

Viasys AVEA

BiCore™ technology integrated into the AVEA ventilator features work of breathing guided ventilator weaning and transpulmonary pressure guidance to allow optimized PEEP titration.

Volume-based capnography with Advanced CO₂ measurements.

Precision gas delivery system that supplies accurate tidal volumes from 2 mL to 2.5 L that automatically compensates for Heliox gas delivery.

Web ID 1879



Philips Respironics V60

Non-Invasive ventilator with an interactive display.

Web ID 2306

Philips Respironics V200

Invasive and Non-Invasive ventilation for pediatrics and adults.

Web ID 2212



Certified Refurbished

Dräger Babylog 8000+

Sensitive synchronization with gentle but precise support for spontaneous breathing reduces the work of breathing and makes the ventilation process much more comfortable for patients.

Harmonious ventilation for small children and the smallest pre-term babies.

Web ID 1883





Baxter InfusO.R.

- Adjust settings easily at any time during infusion.
- Change dial settings with a single hand for an immediate response.
- Size - 9.2 x 4.5 x 2.0 in
- Weight – 2 lbs (with batteries)
- Modify dose without multiple reprogramming steps.
- Minimize risk of calculation error.
- Ideal for administering drugs requiring multiple dosage adjustments.

Easily change infusion parameters and rates with Smart Labels.

Web ID 73



Medfusion 3500

- Programmable dose limits.
- Rapid occlusion detection; Post-occlusion bolus reduction.
- Connectivity with patient data management systems.
- Easy to use and MRI compatible.
- Up to 10 hours internal battery capacity.
- Accepts all syringe sizes cc - 60cc.
- Mounts easily to an IV pole, an infant's isolet or radiant warmer and rotation capabilities.



Intuitive system designed to enhance safety while maintaining simplicity

Web ID 1750

“Smart” programmable safety limits and recordable events.

DRE SP1500 Plus

- Versatile syringe pump features an expanded infusion rate range of up to 1500 mL/h.
- Intuitive and user-friendly operation.
- Simple entry of parameters via numeric keypad.
- Automatic and manual bolus modes.
- Anti-bolus and anti-siphon functions increase patient safety

Web ID 1659

Safe, accurate and easy-to-use.



Alaris 8000 & 8100

Delivers controlled amounts of medications or other fluids to patients through IV, IA, epidural, and other acceptable routes of administration.

Works with Alaris® PC Infusion Pump System; a modular system that builds on a PC-based platform so you can customize infusion delivery by adding and subtracting modules.

Web ID 2187



Hospira Plum A+

Clinicians can choose from a variety of programming options, such as automated piggyback delivery, automated concurrent delivery, and more.

Web ID 2186



Baxter Sigma Spectrum

- Light-weight and compact pump with bright color screen and lithium-ion battery pack.
- Quick search drug libraries assist is easy set-up of infusion.
- Standard IV administration set compatible.
- EMR Integration and Quick-Search MDL Access

Web ID 2405

Aspirators



DRE DM-660

- High performance 1/5 horsepower pump delivers 0 - 26 inHg (0 - 660 mmHg).
- Float shut-off in the canister lid and an Inline Hydrophobic Filter to prevent back flow of fluid as well as airborne contaminants from entering the pump.

Web ID 2107

Gomco Optivac G180

25-550 mm Hg vacuum range, coupled with greater than 30 lpm flow.

Compact size with rechargeable battery that runs for approximately 3 hours at maximum vacuum.

Web ID 1437



DRE DM-530

- High-performance, battery-operated suction unit with battery backup; 1.5 hours of run time.
- Inline Hydrophobic Filter to prevent back flow of fluid as well as airborne contaminants from entering the pump.

Web ID 2392

DRE Comfort Series

Blanket and fluid warmers.

Rapid warming time and uniform heating.

Digital controls and display.

Made of 100% USA Stainless Steel.

Stainless steel or full view glass door, with adjustable shelving.

Two-year warranty included.

9 Different Styles & Sizes
Floor & Countertop Options



Comfort Dual View

Web ID 482



Comfort VC21

Web ID 481



Comfort VC9

Web ID 478



Comfort Dual

Web ID 483



Comfort C5

Web ID 477



Comfort VC5

Web ID 476



Warm Touch 6000

Helps minimize post-surgical complications and improve recovery time.

Intuitive LCD screen and menu function provide all information at a glance.

Web ID 2336

Bair Hugger 505

Snap-fit hose swivels at three points for blanket attachment and positioning.

Web ID 1816



Certified Refurbished

DRE Premio E250

Back up battery.

Side rail controls.

Safe and firm construction from steel with powder coating and platform of ABS plastic with ventilation holes.

Five functions; including CPR function with reset function.

Locking controls prevent unexpected bed movements.

Web ID 2289



The perfect solution for hospitals, surgery centers and teaching facilities.



Hill-Rom Affinity II & III

Web ID 1859




Hill-Rom Total Care

Web ID 1862



Certified Refurbished

Autoclaves



Wide variety of New Tuttnauer Units

Web ID 1834



Factory refurbished Tuttnauer units

Web ID 525



SciCan Statim 5000 G4

Web ID 2355

Stainless Steel


Surgical Scrub Sinks

Single and Double Station Sinks available.

Made with 100% USA stainless steel.



Call for custom sizes at affordable pricing!



I.V. Poles



Mayo Stands



Kick Buckets



Back (Work) Tables

Next Generation Surgical Tourniquet System

AT4 Electronic Tourniquet

Innovative, dual channel electronic tourniquet.

Operates from its own integral air supply and features on-board cuff and hose storage.

Routine automatic diagnostic self-checks at start-up and at regular intervals during use.

Precise regulators provide accurate pressure control at the touch of a button.

Unique IVRA (Biers Block) procedural guidance safety mechanism

Web ID 2401



- Industry leading repair warranties
- Certified component level engineers
- Fast turnaround times
- Tremendous cost savings versus new purchase

Patient Monitors

When it comes to monitors, you want specific answers with real solutions. We can provide a diagnosis and have your equipment reconditioned in just a few easy steps.

Common issues:

- No power
- Parameters failing
- Corrupt hardware
- Faulty NIBP pumps
- Software compatibility
- LCD issues and replacements

Ultrasound Machines & Probes

DRE is your complete source for ultrasound repair, including ultrasound machines and probes. We offer repair services for all major makes and models of ultrasound equipment.

Services Include:

- Fast probe repair and Hotswap programs
- Portable ultrasound depot repair
- On-site engineer scheduling for larger machines
- ISO Certified repair facilities

Infusion and Syringe Pumps

We offer industry leading, ISO certified repair solutions from a wide variety of OEM manufacturers including Abbott, Alaris, Baxter, Braun, Hospira, Medfusion, Sigma and more.

Services Include:

- Extending 'End of Life'
- PM Repairs and Calibrations
- Corrupt Software
- Component level repairs
- Motor and Drive replacements
- 12 month warranty

Defibrillators

DRE offers certified repair on defibrillators with an industry leading warranty. By sending your equipment to DRE you could save significantly compared with a new purchase. We repair most popular models from Philips, Physio Control / Medtronic and Zoll.

Re-Stock your Medical Supplies at

Capnography & CO₂ Nasal Cannulas

Yearly cannula expenses for your monitors can add up quickly.

Call DRE for significant savings on high quality disposable cannulas.



Cautery



Electrodes, Forceps, Grounding Pads, and more.

Monitoring



ECG Leads/Cables, SpO₂ Sensors/Probes, Batteries, Printer Paper, and more.

Surgical Instruments

Available individually and in pre-packed Instrument Sets



Laryngoscopes



Blades, Handles, Lamps, and Complete Sets/Kits.

Shop online for **30,000+ Products** from the Brands You Trust. Masks, Airways, Oxygen Cells, Table Pads, Bougies and more.

DRE Medical Gas Systems



Superior outlets

- Double O-ring seal that doesn't deform and leak when weight is added to the outlet front
- Outlet is capable of supporting 10 lbs. at 2" without leaking
- Outlets incorporate protective vinyl labels that are durable and scratch resistant



Gas delivery management

- Headwalls, consoles and rail systems
- Gas outlets, manifolds and air compressors, alarm systems and pressure switches
- Fittings and hose assemblies
- Emergency oxygen inlet box
- Nitrogen control panel and conversion kits
- Flowmeters and accessories
- Check, relief and ball valves, and zone valve boxes



Medical Zone Valve Boxes

Incorporate dual ported valves. May be piped from left to right or right to left. The installation and pipeline design have been simplified, which can save time and money for the contractor and facility.



Auto Switch Manifold

Ideal for applications requiring automatic switching of cylinders or banks of cylinders.




ANESTHESIA SYSTEMS



RESPIRATORY VENTILATORS



PATIENT TRANSPORT



SURGERY TABLES



SURGERY LIGHTS

European design and quality engineering.

Compare with industry leading brands!

Great for office based procedures:

- Plastic Surgery & Dermatology • Oral Maxillofacial & ENT • Ophthalmology • OB/GYN Specific Models • Diagnostic & General Medicine



DRE Milano E20 Procedure Chair
(p. 18)



1800 Williamson Court
Louisville, KY 40223

PRSR STD
US POSTAGE
PAID
LOUISVILLE, KY
PERMIT #1810

